# visomat

## handy soft

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Wrist Blood Pressure Monitor

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Instructions for use

**UEBE**  
Germany est.1890

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## Table of contents

|          |                                                  |    |
|----------|--------------------------------------------------|----|
| <b>A</b> | <b>Intended use</b>                              | 36 |
| <b>B</b> | <b>Safety instructions</b>                       |    |
| 1.       | Important patient information                    | 37 |
| 2.       | Important technical information                  | 38 |
| <b>C</b> | <b>Operating the unit</b>                        |    |
| 1.       | Unit description                                 | 41 |
| 2.       | Displays                                         | 42 |
| 3.       | Initial start-up of the device                   | 43 |
| 4.       | Set date/time                                    | 43 |
| 5.       | Charge battery                                   | 44 |
| 6.       | Select user                                      | 45 |
| 7.       | Important instructions for use                   | 46 |
| 8.       | Attaching the cuff                               | 47 |
| 9.       | Body posture when taking measurements            | 47 |
| 10.      | Measuring blood pressure                         | 48 |
| 11.      | Classification of the measured values (WHO)      | 48 |
| 12.      | Display of irregular pulse                       | 49 |
| 13.      | Access stored measured values                    | 50 |
| 14.      | Delete stored measured values                    | 51 |
| <b>D</b> | <b>What you should know about blood pressure</b> |    |
| 1.       | Systolic and diastolic blood pressure values     | 52 |
| 2.       | Why you measure different values                 | 52 |
| 3.       | Why should you measure blood pressure regularly? | 52 |
| <b>E</b> | <b>Technical information</b>                     |    |
| 1.       | Failure messages                                 | 53 |
| 2.       | Reset to factory settings                        | 53 |
| 3.       | Customer service                                 | 54 |
| 4.       | Original spare parts and accessories             | 54 |
| 5.       | Technical Inspection / calibration check         | 54 |
| 6.       | Technical data                                   | 55 |
| 7.       | Explanation of symbols                           | 56 |
| 8.       | Disposal                                         | 56 |
| 9.       | Applied standards                                | 57 |
| 10.      | Electromagnetic compatibility (EMC)              | 58 |

## Table of contents

|          |                             |    |
|----------|-----------------------------|----|
| <b>F</b> | <b>Maintaining the unit</b> | 62 |
| <b>G</b> | <b>Warranty</b>             | 63 |

Thank you for choosing the visomat handy soft wrist blood pressure monitor.

These instructions are intended to help the user to use the unit safely and efficiently, and must be kept with the product and forwarded, if applicable.

It is important that you read all the instructions carefully before using the device. Please pay particular attention to the chapter "Important usage information" on page 46.

If you have any remaining questions regarding handling, maintaining, or applying the device after reading the user manual, please contact our visomat customer service (see page 54).

### Content of the package:

Blood pressure monitor visomat handy soft  
 USB cable  
 Charger type L1  
 User manual  
 Blood pressure pass  
 Storage bag

visomat handy soft is an oscillometric blood pressure measurement device for measuring blood pressure and pulse rate at the wrist. The device is intended for adults with a wrist circumference of 13.5 to 21.5 cm. It may only be used indoors.

Measurements are already performed during inflation. The data memory stores the last 60 measurements for both users independently.

The device is powered by an integrated lithium-ion battery.

The device uses the oscillometric method for measuring the blood pressure and pulse rate. After the device has completed its calibration against the ambient air pressure, the automatically controlled inflating process begins. A sensor measures the pressure oscillation that is caused by the expansion and contraction of the arteries at every heart beat. By measuring the amplitude of these pressure waves, the device calculates the systolic and diastolic blood pressure, as well as the pulse rate.

## 1. Important patient information

- The device may only be used at the wrist and not on any other extremities.
- The device may only be used on adults. Measuring the blood pressure of children requires special skills! Consult your doctor in the event that you would like to measure the blood pressure of a child. Under no circumstances should the device be used on an infant.
- The device must not be used on pregnant women or women who suspect that they may be pregnant, as it has not been tested for this patient group. In cases of pre-eclampsia, the device may only be used after consultation with a doctor.
- Do not under any circumstances place the cuff over a critical area, e.g. a wound, aneurysm, etc. or on an arm with an arteriovenous shunt. Risk of injury! Any supply via an intravascular access point (infusion) or other medical monitoring devices could possibly be interrupted.
- In case you are taking medication, ask your doctor about which is the most appropriate time for taking your blood pressure.
- Do not take a blood pressure measurement on the same side as a mastectomy with simultaneous removal of the axilla lymph nodes.
- Discuss the measured values with your doctor. Do not assess the measurement results yourself. Never change the dose of medicines prescribed by your doctor.
- The results of the automatic blood pressure monitor can be distorted by cardiac dysrhythmia, e.g. atrial fibrillation, premature ventricular contraction or arteriosclerosis. Only use the blood measurement device by yourself if you have discussed this with your doctor first.
- Please wait a few minutes between successive measurements, otherwise the blood flow in the arm is interrupted for too long and injuries could occur.

- The device contains small parts that can be swallowed by children. The charging cable entails the risk of suffocation. The device may not be left to children unattended.
- Never touch both the electrically live parts of the device (e.g., charger or charging cable) and the patient at the same time.
- The blood pressure monitor, cuff, and the charger are suitable for a patient environment.
- The device is not suitable for continuously monitoring blood pressure levels during medical emergencies or operations. The blood flow in the arm may be blocked causing numbness, swellings and discolouration of the fingers.
- The device is not suitable for use during ambulance transport.
- The device may not be used simultaneously with high-frequency surgical devices.
- The device is not intended for commercial use.
- The cuff has been tested for skin tolerance according with ISO 10993-5:2009 and ISO 10993-10:2010.
- The device is not suitable for those with an allergy to plastics, especially polyester and nylon.

## 2. Important technical information


- The unit contains sensitive parts and must be protected from extreme temperature fluctuations, humidity, shocks, dust and direct sunlight.
- The device may only be used indoors.
- In case the pressure of the cuff surpasses 300 mmHg, the device will automatically interrupt the inflation process and will vent the cuff. You

may also interrupt the inflation process or measurements at any time by pressing the start/stop key or by removing the cuff.

- If the cuff does not automatically deflate after use due to a defect, immediately undo the hook and loop fastening on the cuff. A cuff pressure above 300 mmHg or a constant pressure above 15 mmHg over a period of more than 3 minutes can cause ecchymotic bleeding.
- The patient is the intended user. The device may only be used for measurements, for charging the battery and for care as described in the user manual. It may not be used for any other purposes than those described in the user manual.
- To avoid inaccurate readings, please adhere to the intended operating and storage conditions. See "Technical data" on page 55.
- The device comes with a safety switch that turns the device off automatically once the voltage surpasses 8 volt or amperage is above 0.5 ampere.
- The charger shields the device from the electric circuit. Pay attention to the fact that the charger can be easily disconnected any time from the electric circuit.
- The charger is part of the medical electrical system. To charge, only use the provided charger type L1 (see accessories on page 54).
- This device may malfunction if used in the vicinity of wireless communication devices such as local wireless networking devices, mobile phones, cordless phones and base stations, radios or other devices with strong electromagnetic fields. Maintain a minimum distance of 30 cm from such devices during use (see "Electromagnetic compatibility" on page 58).
- The device is not category AP/APG. It may not be used in the presence of inflammable anaesthetics with air, oxygen or nitrogen oxide.


- The device comes with a permanently installed rechargeable lithium-ion battery that cannot be changed. Never try to open the device. This may only be done by an authorised expert (see customer service on page 54).
- Never leave the device close to fire or strong heat. This may cause the battery to explode.
- The device may only be used with the original parts and original accessories (see page 54). When using any other parts, the device may be damaged, thereby endangering the user.
- Before use, check that the device is operating reliably and is in proper condition. Do not use the device if it is damaged, as this can lead to injuries or inaccurate measurements.
- If the device malfunctions or is faulty, please refer to the troubleshooting guide starting on page 53, or contact customer service (see page 54).
- In case of unexpectedly hazardous events during usage, contact the visomat customer service immediately (see page 54).

## 1. Unit description

- |                                                 |                                                              |
|-------------------------------------------------|--------------------------------------------------------------|
| 1 Display                                       | 6 MEM key:<br>Access memory<br>Scroll upwards<br>Select user |
| 2 USB Interface                                 |                                                              |
| 3 Cuff (applied part type BF)                   |                                                              |
| 4 Start/Stop button                             | 7 Colour scale of the WHO<br>categorisation                  |
| 5 SET key:<br>Set date/time<br>Scroll downwards |                                                              |

## 2. Displays

- 1 SYS = Systole (top value)
- 2 DIA = Diastole (bottom value)
- 3 PUL 1/min = Pulse
- 4 Current Time
- 5 Date
- 6 User 2
- 7 User 1
- 8 Battery status (see page 44)
- 9 Pulse signal display during measurement / Display of irregular pulse (see page 49)
- 10 WHO categorisation (see page 48)

## 3. Initial start-up of the device

When using the device for the first time, keep the SET key pressed. You will automatically enter the mode for setting the date and time.

If the device cannot be activated, charge the device first and repeat the procedure (see charge battery page 44).

## 4. Set date/time

Before using the device for the first time, set the date and time so that the measurements can be saved with the correct date and time.

Keep the SET key pressed for 3 seconds while the device is shut off. The hour display will start blinking.


By pressing the MEM key, the value increases by 1. Once the correct value has been set, confirm the settings by pressing the SET key.

Now you can enter the time. By pressing the MEM key the value will increase by 1. Confirm the settings by pressing the SET key.

Repeat these steps for setting the day, month, and year.

After completing the date/time settings, the display will show „DONE“ and the device will shut off.

## 5. Charge battery

The device is equipped with a lithium-ion battery that cannot be changed. Upon delivery, the battery is partially charged. Before using the device for the first time, make sure that the device is fully charged.

Only use the provided charger type L1 for charging. Original replacement components and accessories can be found on page 54.

Charge the battery fully

- if the charger icon  blinks on the display during the measurement,
- if the display brightness is decreasing,
- if you turn the device on and the display remains blank.

When the battery symbol starts blinking, you can perform a few more measurements. You should charge the battery as soon as possible.

When the blinking battery symbol  is jointly displayed with Lo (Low), measurements will not be possible. You will first need to charge the battery to perform further measurements.

Connect the charging cable with the USB connector to the blood pressure device and to the charger. Then connect the charger to a suitable socket.

During the charging process, a filling battery icon with a blinking bar will appear on the display. Once the battery is fully charged, the bar will stop blinking.

When the battery is fully discharged, it may take a couple minutes after the charger has been hooked in for the charger display to appear.

After charging, disconnect the charger from the electric circuit and pull the charging cable from the blood pressure monitor.

During the charging process, you cannot measure your blood pressure with the device.

If the battery is fully discharged, charging may take approx. 1-2 hours. When charging a partially discharged battery, the charging time will shorten accordingly.

With a fully charged battery you will be able to perform approx. 75 measurements. The number of measurements can vary slightly.

The battery is permanently installed in the unit and cannot be changed. It is constructed to last approx. 300 charging cycles, which, in the case of an average of 3 measurements per day, amounts to a life cycle of approx. 20 years.

## 6. Select user

The device is designed for 2 users. The results of both users will be saved independently and can be analysed separately, as well.

If you are using the device alone, you may use both users to save your measurements in the morning and the evening separately. In this case, select user 1 for the morning and user 2 for the evening.

Select the user before taking the measurement. With the device switched off, press and hold the MEM button. The user that is currently active is shown.

Change the user by pressing the MEM key until the selected user blinks. Confirm the settings with the SET key.


## 7. Important instructions for use

- To avoid imprecise measurements, please do not conduct measurements:
  - up until 1 hour after any meals or fluid intake,
  - right after consuming tea, coffee, or tobacco,
  - up until 20 minutes after taking a bath,
  - in a cold environment,
  - if you feel increased urinary urgency.
- Measurements must be taken with a resting heartrate. Rest for at least 5 minutes before taking the measurement. Depending on the intensity of any preceding exertion, you may even need to rest for up to an hour.
- Breath deeply 5 to 6 times before starting the measurement.
- It is absolutely necessary to be calm during the measurement. Disturbing movements, vibrations, speaking, and breathing heavily will affect the measurement and can lead to distorted results.
- Pay attention to the display for irregular pulse signals. If necessary, repeat the measurement under more suitable conditions.
- Repeat the measurement after a short break of the blood circuit in the arm, rest for at least 3 minutes for this purpose.
- To attain comparable results, perform the measurement under conditions that are as similar to each other as possible, e.g. always at the same time of day in the same place.
- To exclude the possibility of differences caused by body side it is important to conduct the measurement on one arm only. Please ask your doctor which side is most appropriate for your measurement.
- Blood pressure is not a fixed value. In patients with hypertension, blood pressure may increase or decrease by more than 20 mmHg over a period of a few minutes.

## 8. Attaching the cuff

Clear the wrist from all clothing or watches/jewellery since this may disturb the blood flow in the arm as it may impair the blood pressure at the place of measurement.

Pull the device over your left wrist, palm facing up. The distance between the palm and the cuff's edge must be 1-2 cm.

Now tighten the loose end of the cuff and close it with the Velcro fastener. A cuff that is too loose can cause imprecise measurement values.

## 9. Body posture when taking measurements

Sit at a dining table, lean your back against the backrest of the chair and bend your arm so that your arm is exactly at heart level. Place your feet on the ground and do not cross your legs.

Please pay attention to the following:

- Do not bend your wrist.
- Keep your hand relaxed, not tightened, do not make a fist.
- It is important to keep your arm calm and not to speak during the measurement.

If the cuff is in a higher or lower position than the heart, the measured value might be higher or lower, respectively.

## 10. Measuring blood pressure

After selecting the user and putting on the cuff, press the START/STOP key to start the measurement.

The complete display notification is visible (image 1).

Image 1

After the device has completed its calibration against the ambient air pressure (image 2), the automatically controlled inflating process begins.

Image 2

Measurements are automatically taken during the inflation process (image 3). During measurement, the ♥ symbol blinks to display the heart rate. After the measurements are complete, the results will be shown on the display (image 4).

Image 3

The measured values will be saved automatically.

To switch off the device, press the START/STOP button or wait until the device switches off automatically.

Image 4

## 11. Classification of the measured values (WHO)

The device will classify the measured blood pressure value in accordance with the guidelines of the World Health Organisation (WHO) and the International Society of Hypertension (ISH) of 1999. By looking at the triangle on the display and the colour fields next to display, you will be able to read the classification taken after every measurement.

| Classification                   | Systolic blood pressure<br>= upper value<br>mmHg** | Diastolic blood pressure<br>= lower value<br>mmHg** |
|----------------------------------|----------------------------------------------------|-----------------------------------------------------|
| schwere Hypertonie (Stufe 3)     | ≥ 180*                                             | ≥ 110*                                              |
| mid-level hypertension (level 2) | 160*-179                                           | 100*-109                                            |
| mild hypertension (level 1)      | 140*-159                                           | 90*-99                                              |
| Highly normal                    | 130-139                                            | 85-89                                               |
| Normal                           | 120-129                                            | 80-84                                               |
| Optimal                          | < 120                                              | < 80                                                |

\* It already suffices for one value to be elevated.

\*\* Millimetre hydrargyrum/mercury

Only a doctor can determine your correct blood pressure range and can evaluate if your blood pressure has reached a possibly dangerous level for your case. Discuss your blood pressure values with your doctor.

## 12. Display of irregular pulse

If the device has registered irregular pulses during the measurement, the symbol ♥ for irregular pulse will appear in the display after the measurement. The irregularities may also be caused by cardiac dysrhythmia, disturbing movements, speaking, or by deep breathing.

This symbol will appear if the heart rhythm fluctuates while the device is measuring the systolic and diastolic blood pressure. During the measurement, the device records all pulse intervals and calculates the average value. If two or more pulse intervals deviate from the average by more than 25% or four or more pulse intervals deviate from the average by more than 15%, the ♥ symbol will be displayed to indicate an irregular pulse.

## Operating the unit

If the symbol appears once, it is not usually a reason to worry. If, however, the symbol occurs more frequently, please consult your doctor.

Do not confuse the notification of irregular pulse with the blinking pulse notice which shows the heart rate at every measurement (see page 48).

Please consider that this device is no substitute for a professional examination. It can simply help to detect pulse irregularities at an early stage.

### 13. Access stored measured values

The device holds two measurement storages (one per user) with 60 memory slots each. Measurement results will be automatically stored in the memory. Once 60 values have been saved, the oldest value (no. 60) will be deleted to make room for the newest (no. 1).


Select the user for which you want to access the measured values (see page 45).

In order to access the data, press the MEM key when the device is switched off.

The average value of the last three results for the selected user is indicated with an „A“ (average).

The preceding values can be accessed by pressing the MEM key again (storage number 01). The measured values will alternate between storage number and date/time.

By pressing the MEM and SET key you can scroll through the ascending or descending measured values.


If there are no values stored for a user, the display will show the following:

## Operating the unit

To switch off the device, press the START/STOP button or wait until the device switches off automatically.

### 14. Delete stored measured values

You can delete all values of a user. It is not possible to delete values individually. If you delete the values of one user, the values of the other user will be maintained.

Select the user whose data you wish to delete (see page 45).

Press the MEM key in the shut-off state to get to the memory of the user.

Press the MEM and the SET key simultaneously and keep both pressed for 3 seconds.

The display will show „dEL dONE“ (delete done). The device will shut off automatically afterwards.

## D What you should know about blood pressure

### 1. Systolic and diastolic blood pressure values

The cardiovascular system has the important function of supplying all organs and tissues in the body with sufficient amounts of blood and of transporting metabolites. To do so, the heart contracts and expands at a regular rate of about 60 to 80 times per minute. The pressure of the flowing blood on the artery walls caused by the heart contracting is termed systolic. The pressure in the ensuing relaxation phase, when the heart re-fills with blood, is termed diastolic. During daily measurement you determine both values.

### 2. Why you measure different values

Our blood pressure responds to internal and external influences like a sensitive measuring instrument. Affected by mental, physical and environmental influences, it varies all the time and never remains constant. Reasons for fluctuating blood pressure values can be: moving, speaking, eating, consuming alcohol or nicotine, nervousness, inner tension, emotions, room temperature, recent urination or bowel movements, environmental influences such as movements and sounds, etc. Even changes in the weather and climate can have an effect on your blood pressure. This explains why values measured at the doctor are often higher than those you obtain at home in your usual environment.

### 3. Why should you measure blood pressure regularly?

Even the time of day has an effect on your blood pressure. The values are at their highest during the day. In the course of the afternoon and in the evening, they drop slightly. They are low while you are sleeping, but rise again relatively quickly once you get up.

One-off and irregular measurements therefore say little about your actual blood pressure. A reliable assessment is only possible when measurements are taken regularly. Discuss the measured values with your doctor.

## E Technical information

### 1. Failure messages

| Failure encountered                                                                                                                | Possible cause                                                                | Corrective action                                                                                                           |
|------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| After connecting to the charger, the display will not show the charge indicator.                                                   | Deep discharge of the battery                                                 | Continue charging the device. The charge indicator will appear once the voltage in the battery is high enough. See page 44. |
| Battery notification  + Lo appear in the display. | Battery weak or empty.                                                        | Charge battery                                                                                                              |
| Error message E 1                                                                                                                  | The wristband has not been fastened correctly.                                | Adjust the cuff so it is fastened correctly and repeat the measurement.                                                     |
| Error message E 2                                                                                                                  | The cuff has been fastened too tightly.                                       | Adjust the cuff so it is fastened correctly and repeat the measurement.                                                     |
| Error message E 3                                                                                                                  | The pressure in the cuff is too high.                                         | Relax for a moment and repeat the measurement.                                                                              |
| Error message E 10 or E 11                                                                                                         | Movement or talking detected during the measurement or the pulse is too weak. | Movements can influence your measurement. Relax for a moment and repeat the measurement.                                    |
| Error message E 20                                                                                                                 | No pulse signal detected                                                      | Remove any clothing from your arm and repeat the measurement.                                                               |
| Error message E 21                                                                                                                 | Error during measurement                                                      | Relax for a moment and repeat the measurement.                                                                              |
| Error message EE xx                                                                                                                | Calibration error                                                             | Repeat the measurement. If the error reoccurs, please consult visomat customer service.                                     |
| Error message out                                                                                                                  | Values outside of the measuring range                                         | Relax for a while. Refasten the cuff and repeat the measurement. If the problem persists, consult your doctor.              |

### 2. Reset to factory settings

To reset the device to factory settings, shut the device off and press the MEM and the Start/Stop simultaneously until SHIP appears in the lower edge of the display. Then press the SET key, keeping it pressed until three bars appear on the display. Now the monitor is back to its state as-delivered. To reactivate the monitor, please carry out the procedure for initial start-up (see page 43).

### 3. Customer service

The unit may only be repaired by the manufacturer or by an expressly authorised organisation. Please contact:

UEBE Medical GmbH info@uebe.com  
 Bürgermeister-Kuhn-Straße 22 www.uebe.com  
 97900 Kilsheim, Germany

### 4. Original spare parts and accessories

The following original spare parts and accessories are available from specialist dealers:

- Cuff type H6 (for wrist circumference 13.5-21.5 cm)  
Part. no. 2106001, PZN 11145115
- Charger type L1 with USB charging cable  
Part. no. 2106020, PZN 11145121

Subject to technical modifications.

### 5. Technical inspection / calibration check

Basically, a metrological inspection is recommended at intervals of 2 years. However, professional users in Germany are obligated to comply with the aforementioned according to „Regulation for Operators of Medical Devices“.

This can be implemented either by UEBE Medical GmbH, an authority responsible for metrology, or authorised maintenance services. For this, please observe your national provisions.

Upon request, responsible authorities or authorised maintenance services receive a „Test instruction for metrological inspection“ from the manufacturer.

Please include the manual when delivering the monitor for technical inspection.

Important: No modifications, e.g. opening the unit may be made to this unit without the manufacturer's prior permission.

### 6. Technical data

|                                                                                   |                                                                                                                                  |
|-----------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------|
| Power supply:                                                                     | Built-in lithium-ion battery, 3.7 V, 420 mAh, 5 VDC, 1 A, service life: 300 charge cycles                                        |
| Charger:                                                                          | Type-L1 USB charger, model BLJ06L050100U-V<br>Input 100–240 VAC, 50–60 Hz, max. 0.2 A, output 5 VDC, 1 A                         |
| Display:                                                                          | Digital LCD display, visible field 46.5 mm x 36.5 mm                                                                             |
| Measurement procedure:                                                            | Oscillometric measurement                                                                                                        |
| Reference procedure of clinical testing:                                          | Auscultatory measurement                                                                                                         |
| Pressure display range:                                                           | 0 – 299 mmHg                                                                                                                     |
| Measurement range:                                                                | Systolic: 60 – 230 mmHg<br>Diastolic: 40 – 130 mmHg<br>Pulse: 40 – 199 pulses/min                                                |
| Measurement precision at operation conditions:                                    | Pressure measurement: ± 3 mmHg<br>Pulse measurement: ± 5 %                                                                       |
| Operating conditions:                                                             | Ambient temperature 5 to 40 °C, rel. air humidity 15 to 90 %, air pressure 700 to 1060 hPa                                       |
| Storage and transport conditions:                                                 | Ambient temperature -20 to 60 °C, rel. air humidity up to 93 %, not condensing                                                   |
| Time to adjust from minimum/maximum storage temperature to operating temperature: | 30 minutes                                                                                                                       |
| Cuff:                                                                             | For wrist circumference 13.5-21.5 cm                                                                                             |
| Weight:                                                                           | approx. 103 g                                                                                                                    |
| Dimensions:                                                                       | approx. 80 mm x 73 mm x 13 mm                                                                                                    |
| Accessories:                                                                      | Charger, USB charging cable, user manual, blood pressure diary, storage bag                                                      |
| Mode of operation:                                                                | Continuous operation                                                                                                             |
| Type of the application part:                                                     | BF                                                                                                                               |
| Maximum temperature of the application part:                                      | 42.5°C at an ambient temperature of 40°C                                                                                         |
| Classification:                                                                   | Internal power supply<br>Charger: protection class II                                                                            |
| IP Classification:                                                                | IP 22: Protection against solid foreign bodies with a diameter of more than 12.5 mm; protection against small amounts of liquid. |

|                                |             |
|--------------------------------|-------------|
| Minimum expected service life: | 2 years     |
| Software:                      | Version A01 |

## 7. Explanation of symbols

This product complies with the Council Directive 93/42/EC from 5 September 2007 regarding medical devices and bears the mark CE 0123 (TÜV SÜD Product Service GmbH).

Degree of protection against electric shock: TYPE BF

Please observe instructions for use

Please observe instructions for use

Green Dot

Manufacturer

Reference number = item number

Serial number

Alternating current

Direct current

Equipment class II (double-insulated)

For indoor use only

## 8. Disposal

Batteries and technical appliances must not be disposed of with domestic waste, but should be handed in at the appropriate collection and disposal points.

Please consult visomat customer service if you have questions on disposing the integrated lithium-ion batteries (see page 54).

## 9. Applied standards

|                                 |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|---------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Risk management                 | EN ISO 14971:2012 / ISO 14971:2007 Medical devices - Application of risk management to medical devices                                                                                                                                                                                                                                                                                                                                                                          |
| Labeling                        | EN ISO 15223-1:2016 / ISO 15223-1:2016 Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements                                                                                                                                                                                                                                                                                                 |
| User manual                     | EN 1041:2008 Information supplied by the manufacturer of medical devices                                                                                                                                                                                                                                                                                                                                                                                                        |
| General Requirements for Safety | EN 60601-1:2006+A1:2013 / IEC 60601-1:2005+A1:2012 Medical electrical equipment - Part 1: General requirements for basic safety and essential performance<br>EN 60601-1-11:2015 / IEC 60601-1-11:2015 Medical electrical equipment - Part 1-11: General requirements for basic safety and essential performance - Collateral standard: Requirements for medical electrical equipment and medical electrical systems used in the home healthcare environment                     |
| Electromagnetic compatibility   | EN 60601-1-2:2015 / IEC 60601-1-2:2014 Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances - Requirements and tests                                                                                                                                                                                                                                                    |
| Performance requirements        | EN ISO 81060-1:2012 Non-invasive sphygmomanometers - Part 1: Requirements and test methods for non-automated measurement type<br>EN 1060-3:1997+A2:2009 Non-invasive sphygmomanometers - Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems<br>IEC 80601-2-30:2009+A1:2013 Medical electrical equipment - Part 2-30: Particular requirements for the basic safety and essential performance of automated non-invasive sphygmomanometers |
| Clinical investigation          | EN 1060-4:2004 Non-invasive sphygmomanometers - Part 4: Test procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers<br>ISO 81060-2:2013 Non-invasive sphygmomanometers - Part 2: Clinical validation of automated measurement type                                                                                                                                                                                                     |
| Usability                       | EN 60601-1-6:2010+A1:2015 / IEC 60601-1-6:2010+A1:2013 Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability<br>IEC 62366-1:2015 Medical devices - Part 1: Application of usability engineering to medical devices                                                                                                                                                                          |
| Software life-cycle processes   | EN 62304:2006/AC: 2008 / IEC 62304: 2006+A1:2015 Medical device software - Software life-cycle processes                                                                                                                                                                                                                                                                                                                                                                        |

|                   |                                                                                                                             |
|-------------------|-----------------------------------------------------------------------------------------------------------------------------|
| Bio-compatibility | ISO 10993-1:2009 Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process |
|                   | ISO 10993-5:2009 Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity                         |
|                   | ISO 10993-10:2010 Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization           |

## 10. Electromagnetic compatibility (EMC)

The unit is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use.

Portable and mobile HF communications equipment can affect the unit. Use of the unit in conjunction with non-approved accessories can affect the unit negatively and alter the electromagnetic compatibility. The unit should not be used directly adjacent to or between other electrical equipment.

The unit satisfies the EMC requirements of the international standard IEC60601-1-2. The requirements are satisfied under the conditions described in the tables below.

**Table 1 Guidance and manufacturer's declaration - electromagnetic emissions**

| The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should ensure that it is used in such an environment. |            |                                                                                                                                                                                                                         |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Emissions test                                                                                                                                                                  | Compliance | Electromagnetic environment - guidance                                                                                                                                                                                  |
| RF emissions CIS-PR 11                                                                                                                                                          | Group 1    | The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.                                         |
| RF emissions CIS-PR 11                                                                                                                                                          | Class B    | The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| Harmonic emissions IEC 61000-3-2                                                                                                                                                | Class A    |                                                                                                                                                                                                                         |
| Voltage fluctuations/Flicker emissions IEC 61000-3-3                                                                                                                            | Complies   |                                                                                                                                                                                                                         |

**Table 2 Guidance and manufacturer's declaration - electromagnetic immunity**

| The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment. |                                                                                                                                                                                                                          |                                                                                                                                                                                                                          |                                                                                                                                                |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunity test                                                                                                                                                                   | IEC 60601 test level                                                                                                                                                                                                     | Compliance level                                                                                                                                                                                                         | Electromagnetic environment - guidance                                                                                                         |
| Electrostatic discharge (ESD) IEC 61000-4-2                                                                                                                                     | ±8 kV contact<br>±15 kV air                                                                                                                                                                                              | ±8 kV contact<br>±15 kV air                                                                                                                                                                                              | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %. |
| Electrical fast transient/burst IEC 61000-4-4                                                                                                                                   | ±2 kV for power supply lines<br>±1 kV for input/output lines                                                                                                                                                             | ±2 kV for power supply lines                                                                                                                                                                                             | Mains power quality should be that of a typical commercial or hospital environment.                                                            |
| Surge IEC 61000-4-5                                                                                                                                                             | line(s) to line(s):<br>±1 kV<br>line(s) to earth: ±2 kV<br>100 kHz repetition frequency                                                                                                                                  | line(s) to line(s):<br>±1 kV<br>100 kHz repetition frequency                                                                                                                                                             | Mains power quality should be that of a typical commercial or hospital environment.                                                            |
| Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11                                                                             | 0 % U <sub>i</sub> for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°<br>0 % U <sub>i</sub> for 1 cycle and 70 % U <sub>i</sub> for 25/30 cycles.<br>Single phase at 0°<br>0 % U <sub>i</sub> for 300 cycles | 0 % U <sub>i</sub> for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°<br>0 % U <sub>i</sub> for 1 cycle and 70 % U <sub>i</sub> for 25/30 cycles.<br>Single phase at 0°<br>0 % U <sub>i</sub> for 300 cycles | Mains power quality should be that of a typical commercial or hospital environment.                                                            |
| Power frequency (50/60 Hz) magnetic field IEC 61000-4-8                                                                                                                         | 30 A/m<br>50 Hz/60 Hz                                                                                                                                                                                                    | 30 A/m<br>50 Hz/60 Hz                                                                                                                                                                                                    | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.      |
| NOTE: U <sub>i</sub> is the a.c. mains voltage prior to application of the test level.                                                                                          |                                                                                                                                                                                                                          |                                                                                                                                                                                                                          |                                                                                                                                                |

**Table 3 Guidance and manufacturer's declaration - electromagnetic immunity**

| The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment. |                                                                                                             |                                                                                                             |                                                                                                                                                                                                                                                                                                                   |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Immunity test                                                                                                                                                                   | IEC 60601 test level                                                                                        | Compliance level                                                                                            | Electromagnetic environment - guidance                                                                                                                                                                                                                                                                            |
| Conducted RF IEC 61000-4-6                                                                                                                                                      | 3 V <sub>m</sub> 150 kHz to 80 MHz<br>6 V <sub>m</sub> (in ISM and amateur radio bands)<br>80 % AM at 1 kHz | 3 V <sub>m</sub> 150 kHz to 80 MHz<br>6 V <sub>m</sub> (in ISM and amateur radio bands)<br>80 % AM at 1 kHz | Portable and mobile RF radio equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.<br>Recommended separation distance:<br>$d = 0.35\sqrt{P}$<br>$d = 1.2\sqrt{P}$ |

|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |                          |                          |                                                                                                                                                                                                                                                                                                                                               |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------|--------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Radiated RF IEC 61000-4-3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | 10 V/m, 80 % AM at 1 kHz | 10 V/m, 80 % AM at 1 kHz | d=1.2√P 80 MHz to 800 MHz<br>d=2.3√P 800 MHz to 2.7 GHz<br>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m) b.                                                                                                 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |                          |                          | Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey a, should be less than the compliance level in each frequency range b.<br>Interference may occur in the vicinity of equipment marked with the following symbol:  |
| NOTE1: At 80 MHz and 800 MHz, the higher frequency range applies.<br>NOTE2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |                          |                          |                                                                                                                                                                                                                                                                                                                                               |
| a: Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio stations and TV stations cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.<br>b: Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m. |                          |                          |                                                                                                                                                                                                                                                                                                                                               |

**Table 4 Recommended separation distances between portable and mobile RF communications equipment and the device.**

| The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.                                                                                                                                                |                                                             |                                |                                 |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------|--------------------------------|---------------------------------|
| Rated maximum output power of transmitter W                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 | Separation distance according to frequency of transmitter m |                                |                                 |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             | 150 kHz to 80 MHz<br>d = 3.5√P                              | 80 MHz to 800 MHz<br>d = 1.2√P | 800 MHz to 2.7 GHz<br>d = 2.3√P |
| 0.01                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | 0.12                                                        | 0.12                           | 0.23                            |
| 0.1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | 0.38                                                        | 0.38                           | 0.73                            |
| 1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           | 1.2                                                         | 1.2                            | 2.3                             |
| 10                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          | 3.8                                                         | 3.8                            | 7.3                             |
| 100                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | 12                                                          | 12                             | 23                              |
| For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.<br>NOTE1: At 80 MHz and 800 MHz, the higher frequency range applies.<br>NOTE2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people. |                                                             |                                |                                 |

**Table 5 Guidance and manufacturer's declaration - electromagnetic immunity**

| The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device should assure that it is used in such an environment.                                                                                                                                                                                                                                                                                                                    |              |                                                                 |                                         |                  |            |                         |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------|-----------------------------------------------------------------|-----------------------------------------|------------------|------------|-------------------------|
| Radiated RF IEC 61000-4-3 – Test specifications for enclosure port immunity to RF wireless communications equipment                                                                                                                                                                                                                                                                                                                                                                                                                |              |                                                                 |                                         |                  |            |                         |
| Test frequency MHz                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 | Band a MHz   | Service a                                                       | Modulation b                            | Maximum output W | Distance m | Immunity test level V/m |
| 385                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                | 380 to 390   | TETRA 400                                                       | Pulse modulation b 18 Hz                | 1.8              | 0.3        | 27                      |
| 450                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                | 430 to 470   | GMRS 460, FRS 460                                               | FM c<br>± 5 kHz deviation<br>1 kHz sine | 2                | 0.3        | 28                      |
| 710<br>745<br>780                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | 704 to 787   | LTE Band 13, 17                                                 | Pulse modulation b 217 Hz               | 0.2              | 0.3        | 9                       |
| 810<br>870<br>930                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | 800 to 960   | GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5          | Pulse modulation b 18 Hz                | 2                | 0.3        | 28                      |
| 1720<br>1845<br>1970                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | 1700 to 1990 | GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS | Pulse modulation b 217 Hz               | 2                | 0.3        | 28                      |
| 2450                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | 2400 to 2570 | Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7             | Puls modulation b 217 Hz                | 2                | 0.3        | 28                      |
| 5240<br>5500<br>5785                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               | 5100 to 5800 | WLAN 802.11 a/n                                                 | Pulse modulation b 217 Hz               | 0.2              | 0.3        | 9                       |
| NOTE: If necessary to achieve the immunity test level, the distance between the transmitting antenna and the ME equipment or ME system may be reduced to 1 m. The 1 m distance is permitted by IEC 61000-4-3.<br>a For some services, only the uplink frequencies are included.<br>b The carrier shall be modulated using a 50 % duty cycle square wave signal.<br>c As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case. |              |                                                                 |                                         |                  |            |                         |
| The manufacturer should consider reducing the minimum separation distance, based on risk management and using higher immunity test levels that are appropriate for the reduced minimum separation distance. Minimum separation distance for higher immunity test levels shall be calculated using the following equation:<br>$E = \frac{6}{d} \sqrt{P}$ where P is the maximum output power in W, d is the minimum separation distance in m and E is the immunity test level in V/m.                                               |              |                                                                 |                                         |                  |            |                         |

To avoid imprecise damages or interferences, please consider the following instructions:

- Clean the appliance before and after use with a soft, damp cloth. Do not use abrasive cleaners, benzene, thinners or other strong solvents.
- Avoid dipping the device into water. If the monitor has been dipped into water, clean it with a dry towel.
- Avoid shocks and jolts.
- Avoid washing the cuff. Do not clean cuff or device in a washing machine or dishwasher.
- Never clean the device while it is in use or connected to the charger. Switch off the device and disconnect it from the charger before cleaning.
- Store the device in a dry, cool, well-ventilated place that is protected from direct sunlight.
- Protect the device from strong temperature fluctuations and dusty environments.
- Take care when storing the device. Ensure that no heavy objects are resting on the device or the cuff.
- If the unit has been stored at temperatures below 10 °C or above 40°C, let it stand at room temperature for at least 30 minutes before use.

The device has been manufactured and tested with great care. However, in the unlikely event of a defect being detected after delivery, we provide warranty in accordance with the following terms and conditions:

During the warranty period of 5 years from the date of purchase we reserve the right either to repair any such defect at our expense or to supply a perfect replacement unit.

Excluded from the warranty are parts subject to normal wear and tear as well as damage caused by non-compliance with the instructions for use, improper handling (e.g. unsuitable power sources, breakages, leaking batteries) and/or disassembly of the unit by the purchaser. Furthermore, no claims for damages against us are substantiated by the warranty.

Warranty claims can only be advanced in the warranty period and by presenting proof of purchase. In the event of a warranty claim, the unit must be sent to the following address together with the proof of purchase and a description of the complaint:

UEBE Medical GmbH  
Service-Center  
Zum Schlag 18  
97877 Wertheim, Germany

The cost of returning the unit to our factory shall be borne by the sender. UEBE shall refuse to accept return deliveries that have not been paid for by the sender.

The statutory claims and rights of the buyer against the seller (claims for defect, manufacturer's liability, for example) are not restricted by this warranty.

**Please note: In the event of a warranty claim it is essential to attach the proof of purchase.**

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© Copyright 2020 UEBE Medical GmbH