# visomat

## comfort soft *BT*

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## Table of contents

<b>A</b>	<b>Intended use</b>	34
<b>B</b>	<b>Safety instructions</b>	
	1. Important patient information	35
	2. Important technical information	36
<b>C</b>	<b>Operating the device</b>	
	1. Device description	39
	2. Display	40
	3. Important usage information	41
	4. Inserting/replacing batteries	42
	5. Setting the date/time	42
	6. Bluetooth data transfer	43
	7. Attaching the cuff	43
	8. Body position during the measurement	45
	9. Selecting the user	45
	10. Measuring blood pressure	46
	11. Classification of the measured values (WHO)	46
	12. Display of irregular pulse	47
	13. Using the memory	48
<b>D</b>	<b>What you should know about blood pressure</b>	
	1. The systolic and diastolic blood pressure value	50
	2. Why you measure different values	50
	3. Why should you measure blood pressure regularly?	50
<b>E</b>	<b>Technical information</b>	
	1. Error messages	51
	2. Customer service	51
	3. Technical data	52
	4. Original spare parts and accessories	53
	5. Metrological inspection	53
	6. Explanation of symbols	54
	7. Disposal	55
	8. Electromagnetic compatibility (EMC)	55

<b>F</b>	<b>Maintaining the device</b>	58
<b>G</b>	<b>Warranty</b>	59

Thank you for choosing the visomat comfort soft BT upper arm blood pressure monitor,

These instructions for use are intended to help the user operate the device safely and efficiently, and must be kept with the product and forwarded, if applicable. It is important that you read all the instructions carefully before using the device. Please pay particular attention to the chapter "Important usage information" on page 41.

If you have any remaining questions regarding handling, maintaining, or applying the device after reading the user manual, please contact our visomat customer service (see page 51).

### Content of the package:

Blood pressure monitor visomat comfort soft BT  
Cuff type WRT for arm circumference of 22-42 cm  
Transport and storage pouch  
4 batteries AAA 1.5 V  
Instructions for use

visomat comfort soft BT is an oscillometric blood pressure measurement device for measuring blood pressure and pulse rate at the upper arm. The device is intended for adults with a upper arm circumference of 22 to 42 cm. It may only be used indoors.

Measurements are already performed during inflation. The data memory stores the last 250 measurements for both users independently.

The device uses the oscillometric method for measuring the blood pressure and pulse rate. After the device has completed its calibration against the ambient air pressure, the automatically controlled inflating process begins. A sensor measures the pressure oscillation that is caused by the expansion and contraction of the arteries at every heart beat. By measuring the amplitude of these pressure waves, the device calculates the systolic and diastolic blood pressure, as well as the pulse rate.

## 1. Important patient information

- The device may only be used at the upper arm and not on any other extremities.
- The device may only be used on adults. Measuring the blood pressure of children requires special skills! Consult your doctor in the event that you would like to measure the blood pressure of a child. Under no circumstances should the device be used on an infant.
- The device must not be used on pregnant women or women who suspect that they may be pregnant, as it has not been tested for this patient group. In cases of pre-eclampsia, the device may only be used after consultation with a doctor.
- Do not under any circumstances place the cuff over a critical area, e.g. a wound, aneurysm, etc. or on an arm with an arteriovenous shunt. Risk of injury! Any supply via an intravascular access point (infusion) or other medical monitoring devices could possibly be interrupted.
- In case you are taking medication, ask your doctor about which is the most appropriate time for taking your blood pressure.
- Discuss the measured values with your doctor. Do not assess the measurement results yourself. Never change the dose of medicines prescribed by your doctor.
- Do not take a blood pressure measurement on the same side as a mastectomy with simultaneous removal of the axilla lymph nodes.
- The results of the automatic blood pressure monitor can be distorted by cardiac dysrhythmia, e.g. atrial fibrillation, premature ventricular contraction or arteriosclerosis. Only use the blood measurement device by yourself if you have discussed this with your doctor first.
- Ensure that the air hose is not kinked. A kinked air hose may prevent the cuff from deflating, thus interrupting the blood flow in the arm for too long.

## Safety instructions

- Please wait a few minutes between successive measurements, otherwise the blood flow in the arm is interrupted for too long and injuries could occur.
- The device contains small parts that could be swallowed by children. The air hose poses a strangulation hazard. Do not leave the device unattended with children.
- Never simultaneously touch the device's live electrical parts and the patient.
- The blood pressure monitor and the cuff are suitable for a patient environment.
- The device is not suitable for continuously monitoring blood pressure levels during medical emergencies or operations. The blood flow in the arm may be blocked causing numbness, swellings and discoloration of the fingers.
- The device is not suitable for use during ambulance transport.
- The device may not be used simultaneously with high-frequency surgical devices.
- The device is not intended for commercial use.
- The device is not suitable for those with an allergy to plastics, especially polyester and nylon.

## 2. Important technical information

- The unit contains sensitive parts and must be protected from extreme temperature fluctuations, humidity, shocks, dust and direct sunlight.
- The device may only be used indoors.

## Safety instructions

- If the unit has been stored at temperatures below 10 °C or above 40 °C, let it stand at room temperature for at least 30 minutes before use.
- In case the pressure of the cuff surpasses 300 mmHg, the device will automatically interrupt the inflation process and will vent the cuff. You may also interrupt the inflation process or measurements at any time by pressing the start/stop key or by removing the cuff.
- If the cuff does not automatically deflate after use due to a defect, immediately undo the hook and loop fastening on the cuff. A cuff pressure above 300 mmHg or a constant pressure above 15 mmHg over a period of more than 3 minutes can cause ecchymotic bleeding.
- The patient is the intended user. The device may only be used for measurements, for changing the battery and for care as described in the user manual. It may not be used for any other purposes than those described in the user manual.
- To avoid inaccurate readings, please adhere to the intended operating and storage conditions. See "Technical data" on page 52.
- This device may malfunction if used in the vicinity of wireless communication devices such as local wireless networking devices, mobile phones, cordless phones and base stations, radios or other devices with strong electromagnetic fields. The device itself may also cause interference with other nearby electronic devices. During use, maintain a minimum distance of 1 m from such devices (see "Electromagnetic compatibility" on page 55).
- At the end of the measurement process, hold the device at least 20 cm away from the body while the data is transferred.
- The device is not category AP/APG. It may not be used in the presence of inflammable anaesthetics with air, oxygen or nitrogen oxide.
- Always replace all batteries simultaneously during battery replacement. Do not mix new and old batteries or batteries of different types.


- You will require 4 x 1.5 V batteries. Rechargeable batteries only have a voltage of 1.2 V and are therefore not suitable.
- The mains adapter is part of the medical electrical system. When using the device with a mains adapter, please only use the visomat mains adapter, which has been specially tested for medical devices. (See "Accessories" on page 53).
- The power cord is classed as an isolator used to disconnect the device from the electricity supply. Make sure that the power cord is always easily accessible so that the device can be disconnected from the electricity supply in an emergency.
- If the device is not to be operated for an extended period of time, please remove the batteries. In principle, all batteries can leak.
- Do not throw batteries onto a fire. They may explode or leak.
- The device may only be used with the original parts and original accessories. When using any other parts, the device may be damaged, thereby endangering the user.
- Never open or modify the device (except for replacing the batteries).
- Before use, check that the device is operating reliably and is in proper condition. Do not use the device if it is damaged, as this can lead to injuries or inaccurate measurements.
- The inflating and measuring procedure can be interrupted by pressing the Start/Stop button or by removing the cuff. The device then stops the inflation procedure and deflates the cuff.
- If the device malfunctions or is faulty, please refer to the troubleshooting guide starting on page 51, or contact customer service (see page 51).
- Any serious adverse effects linked to the product must be reported to the manufacturer and the relevant authorities in the member state in which the user/patient resides.

## 1. Device description

- |  |   |
|--|---|
| <b>1</b> Illuminated display                 | <b>6</b> Settings button (S = Set)                            |
| <b>2</b> Cuff (applied part type BF)         | <b>7</b> Start/Stop button                                    |
| <b>3</b> Connecting socket for mains adapter | <b>8</b> Battery compartment (on the underside of the device) |
| <b>4</b> User selection switch               | <b>9</b> Colour scale of the WHO categorisation               |
| <b>5</b> Memory button (M = Memory)          | <b>10</b> Air tube connector                                  |

## 2. Display

- |  |  |
|--|--|
| <b>1</b> Memory display                | <b>8</b> Cuff deflated                                   |
| <b>2</b> SYS = Systole (top value)     | <b>9</b> Pulse signal display/irregular pulses (page 47) |
| <b>3</b> DIA = Diastole (bottom value) | <b>10</b> Battery indicator (page 42)                    |
| <b>4</b> PUL 1/min = Pulse             | <b>11</b> Average measurement display                    |
| <b>5</b> Date/time                     | <b>12</b> User B   |
| <b>6</b> WHO classification (page 46)  | <b>13</b> User A   |
| <b>7</b> Bluetooth data transfer       |  |

## 3. Important usage information

- To avoid imprecise measurements, please do not conduct measurements:
  - up until 1 hour after any meals or fluid intake,
  - right after consuming tea, coffee, or tobacco,
  - up until 20 minutes after taking a bath,
  - in a cold environment,
  - if you feel increased urinary urgency.
- Measurements must be taken with a resting heartrate. Rest for at least 5 minutes before taking the measurement. Depending on the intensity of any preceding exertion, you may even need to rest for up to an hour.
- Breath deeply 5 to 6 times before starting the measurement.
- It is absolutely necessary to be calm during the measurement. Disturbing movements, vibrations, speaking, and breathing heavily will affect the measurement and can lead to distorted results.
- Pay attention to the display for irregular pulse signals.
- Repeat the measurement after a short break of the blood circuit in the arm, rest for at least 3 minutes for this purpose.
- To attain comparable results, perform the measurement under conditions that are as similar to each other as possible, e.g. always at the same time of day in the same place.
- To exclude the possibility of differences caused by body side it is important to conduct the measurement on one arm only. Please ask your doctor which side is most appropriate for your measurement.
- Blood pressure is not a fixed value. In patients with hypertension, blood pressure may increase or decrease by more than 20 mmHg over a period of a few minutes.

#### 4. Inserting/replacing batteries

- Open the battery compartment on the underside of the device.
- Insert new batteries. Check that the polarity is correct (see markings in the battery compartment).
- Close the battery compartment again.

Replace the batteries

- when the display shows the battery symbol  alongside the text "Lo" (Low).
- the display gets dimmer.
- if you switch on the device and nothing appears on the display.

#### 5. Setting the date/time

Before using the device for the first time, set the date and time so that the measurements can be saved with the correct date and time.

With the device switched off, briefly press the S button. The current time is displayed. Press the S button again and hold for 3 seconds. The year starts to flash (Figure 1).

By pressing the M key, the value increases by 1. Once the correct value has been set, confirm the settings by pressing the S key.

Repeat these steps to enter the month and day (Figure 2) and the hours and minutes (Figure 3).

Figure 1

Figure 2

Figure 3

Once you have completed the basic settings, "dONE" appears on the display (Figure 4). The settings that you made now appear again on the display. The device then switches off automatically.

With the device switched off, if you briefly press the S button, the time appears on the display.

Figure 4

#### 6. Bluetooth data transfer

The device has a Bluetooth interface which enables you to transmit the measurements to your smartphone.

If you wish to use this function, please download the "visomat Blood Pressure Diary" app from the App Store or Google Play.

After installation, open the visomat Blood Pressure Diary app and register as a user. To do this, follow the instructions in the app.

To pair the device, please follow the instructions in the app. Pairing can only be performed after a measurement has been taken and while the Bluetooth symbol  is flashing in the display.

If you do not wish to use the data transmission and evaluation function, you can also use the device without a smartphone.

#### 7. Attaching the cuff

Before applying the cuff, please check that the circumference of your upper arm is within the range specified on the cuff. An incorrect cuff size could result in inaccurate readings.

## Operating the device

- Expose your upper arm. Under no circumstances must clothing or watches/jewellery hinder the blood flow in the arm, as this affects the blood pressure at the measuring point.
- Push the cuff onto the arm until the lower edge of the cuff is 2-3 cm above the crook of the arm (Figure 1).


Figure 1

- When placed on the left arm, the air hose runs down the middle of the crook of the arm to the device. This automatically positions the artery marking centrally above the pulse position (Figure 2).


Figure 2

- When placed on the right arm, the cuff must be rotated to the left until the artery marking is on the pulse position. The air hose then runs along the inner side of the upper arm (Figure 3).


Figure 3

- The cuff should not be too tight. You should be able to insert two fingers between the arm and the cuff. Tip: Angle the arm slightly. Slightly tense the upper arm muscle. This increases the circumference of the arm slightly.


Figure 4

- Now pull the free end of the cuff tight and close the Velcro fastener (Figure 4).
- Check that the measurement arrow on the cuff is within the "Marking for arm circumference" on the cuff edge (Figure 5).


Figure 5

- Connect the plug of the cuff to the cuff connection on the left-hand side of the device. Ensure that the plug is fully inserted into the device.

## Operating the device

### 8. Body position during the measurement

- Sit at a table (preferably at the height of a dining table, not a coffee table).
- Sit with your back against the backrest of the chair.
- Rest your entire lower arm on the table with the palm facing upwards.
- Put your feet on the floor and do not cross your legs.

Ensure that the measuring point on the upper arm is approximately level with your heart. If the measuring point is lower than the heart, the measured values may be higher. If the measuring point is higher than the heart, the measured values may be lower.

### 9. Selecting the user

The device is designed for two users. The values for the two users are stored separately and can also be evaluated separately.

You may also use the two user memories to evaluate measurements taken in morning and those taken in the evening separately. To do this, carry out the measurement in the morning using user A and in the evening using user B.

Move the user selection switch to select a user. Switch to the right = User A, switch to the left = User B.

Before taking the measurement, check that the correct user is displayed.

Now carry out the measurement.

## 10. Measuring blood pressure

- Press the START/STOP button to start measuring your blood pressure. The device carries out the measurement fully automatically and saves the readings.
- The complete display notification is visible (image 1).
- After the device has completed its calibration against the ambient air pressure (Figure 2), the automatically controlled inflation process starts.
- Measurements are automatically taken during the inflation process (image 3). During measurement, the ♥ symbol blinks to display the heart rate. After the measurements are complete, the results will be shown on the display (image 4).
- Once the measurement process is complete, the device looks for a nearby Bluetooth connection to transfer the readings. While this is taking place, the Bluetooth symbol ☎ flashes on the display. Once the data has been successfully transferred, the device waits for 5 seconds then switches itself off.
- If you do not use the Bluetooth function, ignore the ☎ symbol. To switch off the device, press the START/STOP button or wait until the device switches off automatically after 3 minutes.


Figure 1

Figure 2

Figure 3

Figure 4

## 11. Classification of the measured values (WHO)

The device will classify the measured blood pressure value in accordance with the guidelines of the World Health Organisation (WHO) and the In-

ternational Society of Hypertension (ISH) of 1999. By looking at the bars on the display and the colour fields next to display, you will be able to read the classification taken after every measurement.

Classification	Systolic pressure = top value mmHg	Diastolic pressure = bottom value mmHg
strong hypertonia (level 3)	≥ 180	≥ 110
mid-level hypertonia (level 2)	160-179	100-109
mild hypertonia (level 1)	140-159	90-99
High-normal	130-139	85-89
Normal	120-129	80-84
Optimal	< 120	< 80

Only a doctor can determine the right blood pressure range for you and assess whether your blood pressure has reached a dangerous level for you. Discuss your blood pressure values with your doctor.

**Never change the dose of medicines prescribed by your doctor!**

## 12. Display of irregular pulse

If, after a measurement, the irregular pulse symbol ♥ appears on the display, the device has recorded pulse anomalies during the reading. This can be triggered by cardiac arrhythmias, movement, talking or even by deep breathing.

This symbol will appear if the heart rhythm fluctuates while the device is measuring the systolic and diastolic blood pressure. During the measurement, the device records all pulse intervals and calculates the average value. If two or more pulse intervals deviate from the average by more than 25%, or four or more pulse intervals deviate from the average by more than 15%, the ♥ symbol will be displayed to indicate an irregular pulse.

If the symbol appears, this is usually no cause for concern. However, if the symbol appears on a regular basis, consult a doctor for advice.

Please note that using this device is not a substitute for a medical examination. It is intended only to help identify pulse irregularities at an early stage.

Do not confuse the irregular pulse display  with the flashing pulse display  that indicates the heart rate during every measurement (see page 46).

### 13. Using the memory

The device stores the measured values in two memories (one per user). Each memory can store up to 250 readings and the average value. The measured results are automatically saved in the memory. If more than 250 measured values are saved, the oldest measured value (no. 250) is deleted so that the latest value (no. 1) can be recorded.

#### Access stored measurement values

Move the user selection switch to choose the user whose data you wish to view. Switch to the right = User A, switch to the left = User B.

With the device switched off, press the M button once.

If 3 or more readings are stored in the memory, the average value of the last 3 results is displayed next to "AVG" (Average) (Figure 1).

Figure 1

Use the M and S buttons to scroll up or down through the individual readings. If you press and hold the M button, the list scrolls quickly in groups of 10 readings.


Figure 2


Figure 3

The readings appear on the display next to the alternating date and time (Figure 2, Figure 3).

To switch off the device, press the START/STOP button or wait until the device switches off automatically after 1 minute.

If there are no values stored, the display will show the following (Figure 4):

Figure 4

#### Delete saved values

You can delete all values of a user. It is not possible to delete values individually. If you delete the values of one user, the values of the other user will be maintained.

Move the user selection switch to choose the user whose data you wish to view. Switch to the right = User A, switch to the left = User B.

Figure 5

With the device switched off, press the memory button M to access the memory.

Press and hold the S button for 3 seconds. "dEL ALL" (Delete All, Figure 5) appears on the display.

Press the S button again to confirm that you want to delete the data. The display will show "dEL dONE" (delete done, figure 6). The device will shut off automatically afterwards.

Figure 6

## D What you should know about blood pressure

### 1. The systolic and diastolic blood pressure value

The circulatory system has the important task of supplying all organs and tissues in the body with sufficient blood, and of removing metabolites. The heart does this by contracting and expanding at a regular rhythm of around 60-80 times per minute. The pressure created by the blood flow on the artery walls when your heart contracts is referred to as systole. The pressure in the subsequent relaxation phase when the heart fills with blood again is referred to as diastole. In your daily measurement, you will measure both values.

### 2. Why you measure different values

Our blood pressure responds to internal and external influences like a sensitive measuring instrument. Affected by mental, physical and environmental influences, it varies all the time and never remains constant. Reasons for fluctuating blood pressure values can be: moving, speaking, eating, consuming alcohol or nicotine, nervousness, inner tension, emotions, room temperature, recent urination or bowel movements, environmental influences such as movements and sounds, etc. Even changes in the weather and climate can have an effect on your blood pressure.

This explains why values measured at the doctor are often higher than those you obtain at home in your usual environment.

### 3. Why should you measure blood pressure regularly?

Even the time of day has an effect on your blood pressure. The values are at their highest during the day. In the course of the afternoon and in the evening, they drop slightly. They are low while you are sleeping, but rise again relatively quickly once you get up.

One-off and irregular measurements therefore say little about your actual blood pressure. A reliable assessment is only possible when measurements are taken regularly. Discuss the measured values with your doctor.

## E Technical information

### 1. Error messages

Error	Possible cause	Remedy
When you turn the device on, the display is blank or dim.	Batteries low or empty	Insert new batteries
	Batteries not installed correctly	Install the batteries correctly. Check the polarity markings in the battery compartment.
	Mains adapter not connected correctly	Connect the mains adapter correctly
The display shows the battery indicator  + Lo	Batteries low or empty	Insert new batteries
Error message E01	Cuff is too tight or too loose	Adjust the cuff so it is fastened correctly and repeat the measurement.
Error message E02	Movement detected during measurement.	Relax for a moment and repeat the measurement.
Error message E03	No pulse signal detected	Loosen any clothing on your arm and repeat the measurement
Error message E04	Error during measurement	Relax for a moment and repeat the measurement.
Error message EE xx	Calibration error	Repeat the measurement. If the error reoccurs, please consult UEBE customer service.
Error message out	Readings are outside the measurement range	Relax for a moment, remove then replace the cuff on your arm and repeat the measurement. If this message appears again, contact your doctor.

### 2. Customer service

Device repairs may only be performed by the manufacturer or an expressly authorised body. Please contact:

UEBE Medical GmbH  
 Bgm.-Kuhn-Str. 22  
 97900 Kulsheim, Germany  
 info@uebe.com, www.uebe.com

### 3. Technical data

Power supply:	On battery power: 6 V DC, 4 x 1.5 V AAA batteries
	Optional power via mains adapter (not included): visomat mains adapter, output 6 V DC, minimum 500 mA
Display:	Digital LCD display, visible field 60.5 mm x 92.5 mm
Measurement method:	Oscillometric measurement
Reference method of the clinical test:	Auscultatory measurement
Pressure display range:	0 – 299 mmHg
Measurement range:	Systolic: 60-230 mmHg Diastolic: 40-130 mmHg Pulse measurement: 40-199 pulses/minute
Measurement precision at operation conditions:	Pressure measurement: $\pm 3$ mmHg Pulse measurement: $\pm 5$ %
Operating conditions:	Ambient temperature 5 °C to 40 °C, relative humidity 15% to 90% non-condensing, air pressure 700 to 1060 hPa
Conditions for storage and transport:	Ambient temperature -20 to 60 °C, rel. air humidity up to 93 %, not condensing
Time to adjust from minimum/maximum storage temperature to operating temperature:	30 minutes
Cuff:	For upper arm circumference 22-42 cm
Weight:	approx. 246 g without batteries and cuff
Size:	approx. 140 mm x 130 mm x 49.7 mm
Operating mode:	Continuous operation
Type of the application part:	BF
IP rating:	IP 21: Protection against solid foreign bodies with a diameter of more than 12.5 mm; protection against small amounts of liquid.
Expected service life:	2 years
Classification:	Internal power supply using batteries
Software:	Version A01

### 4. Original spare parts and accessories

The following original spare parts or accessories are available through specialist retailers:

Universal cuff 22-42 cm, type WRT  
Part no. 2406001, PZN-13424865

visomat mains adapter A2  
Part no. 2400020, PZN-11870684

Subject to technical modifications.

### 5. Metrological inspection

Basically, a metrological inspection is recommended at intervals of 2 years. However, professional users in Germany are obligated to comply with the aforementioned according to "Regulation for Operators of Medical Devices".

This can be performed either by UEBE Medical GmbH, an authority responsible for metrology or an authorised maintenance service. Please refer to your national regulations.

Upon request, responsible authorities or authorised maintenance services receive a "Test instruction for metrological inspection" from the manufacturer.

Please only submit the device for metrological inspection together with the cuff and the instructions for use.

Caution: No modifications such as opening or repairing the device may be made without the manufacturer's permission.

## 6. Explanation of symbols

This product complies with Council Directive 93/42/EEC of 5 September 2007 concerning medical devices and bears the mark CE 0123 (TÜV SÜD Product Service GmbH).

Consult the instructions for use

Caution: Observe safety instructions and warnings

Manufacturer

Medical Device

Reference number = item number

Lot number

Serial number

Unique Device Identification

Temperature limitation

Humidity limitation

Keep dry

Protection against solid foreign bodies with a diameter of more than 12.5 mm; protection against small amounts of liquid.

Protection class of the applied part against electric shock: TYPE BF

Direct current

Protection class II (double insulation)

For indoor use only

Polarity inside plus

## 7. Disposal

Technical equipment and batteries do not belong in household waste. They must be disposed of at appropriate collection and disposal points.

## 8. Electromagnetic compatibility (EMC)

The unit is an electrical medical product and is subject to special precautionary measures with regard to EMC which must be published in the instructions for use.

Portable and mobile HF communications equipment can affect the unit. Use of the unit in conjunction with non-approved accessories can affect the unit negatively and alter the electromagnetic compatibility. The unit should not be used directly adjacent to or between other electrical equipment.

The unit satisfies the EMC requirements of the international standard IEC60601-1-2. The requirements are satisfied under the conditions described in the tables below.

### Guidance and manufacturer's declaration - electromagnetic emissions

Emissions test	Compliance
RF emissions CISPR 11	Group 1
RF emissions CISPR 11	Class B
Harmonic emissions IEC 61000-3-2	Class A
Voltage fluctuations/Flicker emissions IEC 61000-3-3	Complies

## Guidance and manufacturer's declaration - electromagnetic immunity

Immunity test	IEC 60601 test level	Compliance level
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	±8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines 100 kHz repetition frequency	±2 kV for power supply lines N/A 100 kHz repetition frequency
Surge IEC 61000-4-5	±0.5 kV, ±1 kV differential mode ±0.5 kV, ±1 kV, ±2 kV common mode	±0.5 kV, ±1 kV differential mode N/A
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0 % U <sub>n</sub> for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U <sub>n</sub> for 1 cycle and 70 % U <sub>n</sub> for 25/30 cycles. Single phase at 0° 0 % U <sub>n</sub> for 250/300 cycles	0 % U <sub>n</sub> for 0.5 cycle at 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0 % U <sub>n</sub> for 1 cycle and 70 % U <sub>n</sub> for 25/30 cycles. Single phase at 0° 0 % U <sub>n</sub> for 250/300 cycles
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m 50 Hz/60 Hz	30 A/m 50 Hz/60 Hz
Conducted RF IEC 61000-4-6	3 Vm 150 kHz to 80 MHz 6 Vm (in ISM and amateur radio bands) 150 kHz to 80 MHz 80 % AM at 1 kHz	3 Vm 150 kHz to 80 MHz 6 Vm (in ISM and amateur radio bands) 150 kHz to 80 MHz 80 % AM at 1 kHz
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz	10 V/m 80 MHz to 2.7 GHz 80 % AM at 1 kHz

NOTE: U<sub>n</sub> is the a.c. mains voltage prior to application of the test level.

## Guidance and manufacturer's declaration - electromagnetic immunity

Radiated RF IEC 61000-4-3 – Test specifications for enclosure port immunity to RF wireless communications equipment							
Test frequency MHz	Band MHz	Service	Modulation	Maximum output W	Distance m	Immunity test level V/m	Compliance level V/m
385	380 to 390	TETRA 400	Pulse modulation 18 Hz	1.8	0.3	27	27
450	430 to 470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0.3	28	28
710	704 to 787	LTE Band 13, 17	Pulse modulation 217 Hz	0.2	0.3	9	9
745							
780							
810 870	800 to 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0.3	28	28
930							
1720	1700 to 1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation 217 Hz	2	0.3	28	28
1845							
1970							
2450	2400 to 2570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450 LTE Band 7	Puls modulation 217 Hz	2	0.3	28	28
5240 5500 5785	5100 to 5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0.2	0.3	9	9

To avoid imprecise damages or interferences, please consider the following instructions:

- Clean the appliance before and after use with a soft, damp cloth. Do not use abrasive cleaners, benzene, thinners or other strong solvents.
- Do not allow the device come into contact with water. However, if it does come into contact with water, wipe the device with a dry cloth.
- Avoid shocks and jolts.
- Do not try to clean the cuff with water. Never immerse it in water. Do not clean the device in a washing machine or dishwasher.
- Never clean the device while it is in use. Switch off the device before cleaning.
- Store the device in a dry, cool, well-ventilated place that is protected from direct sunlight.
- Protect the device from strong temperature fluctuations and dusty environments.
- Take care when storing the device. Ensure that no heavy objects are resting on the device or the cuff.

The device has been manufactured and tested with all due care. Nevertheless, in the event of defects upon delivery, we provide a warranty under the following conditions:

During the warranty period of 5 years from the date of purchase, we remedy such defects at our discretion and at our expense in our plant through repair work or replacement delivery of a defect-free device.

The warranty does not cover normal wear and tear of wearing parts or damage caused by failure to observe the instructions for use, improper handling (e.g. unsuitable power sources, breakage, leaking batteries) and/or disassembly of the device by the buyer. Furthermore, the warranty does not constitute grounds for asserting claims for damages against us.

Warranty claims can only be asserted during the warranty period and upon presentation of the proof of purchase. In the case of a warranty claim, the device along with the proof of purchase and a description of the complaint must be sent to:

UEBE Medical GmbH  
Service-Center  
Bgm.-Kuhn-Str. 20  
97900 KÜlsheim, Germany

The cost of returning the device to our factory shall be borne by the sender. Complaints that are sent back without prepayment will not be accepted by UEBE.

The statutory claims and rights of the buyer against the seller (for example, defect claims, producer liability) are not restricted by this warranty.

Please note: In the case of a warranty claim, please make sure to enclose the proof of purchase.

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© Copyright 2023 UEBE Medical GmbH